Woods lamp

Product code:1004

User manual of Woods
lamp

Lux Kozmetika d.o.o.

WARNING

1. Power supply must be cut off when changing electrical components or maintaining and cleaning machines.
2. The plug must be plugged into an outlet that is properly installed and grounded, to prevent the risk of fire and electric shock.
3. Do not allow children to play with and place the machine in a place where children can avoid touching.
4. Please do not place the machine in damp places such as bathroom to avoid electric shock or other failure.
5. Please do not place the machine in places where painting/coating or oxygen providing during operation.
6. Unless being monitored or instructed by personnel who can be responsible for their personal safety, the people who have physical obstacles, tactile disorders, mental disorders, lacking of experience, common sense or children, are not allow to use the machine.
7. Please contact distributor, service center or related personal with appropriate skill and experience for repairing if the plug is damaged, in case of any dangerous.
8. Maintenance should be processed by trained professionals.
9. This product can only be used by professional beauty practitioners, please pay close attention to working condition of the machine during operating.

Instruction

Ultraviolet ray lamp is made from general ultraviolet ray which go through the glass-filter with nickel. Wood's lamp was invented by American physical scientist Robert Williams Wood's. The main medical usage is to diagnose the unfavorable condition of skin, such as rash.

Usage

Beautician use wood's lamp to examine the condition of skin. Different material irradiated with dark purple ray of wood's lamp will appear ray with different colors. The wood's lamp is supposed to be used in a dark room. The purple ray can be used to examine the scarf skin of human body by beautician. What's more, the beautician can do some appropriate skin care for the customer. The invisible skin blemish can be displayed by wood lamp. The displayed shape and color of different skin would be different. For example, if the skin is too thick, the displayed ray should be bright white. The following cases are some of the situation under wood's lamp.

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